Fuego^{SCS} basic

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The Fuego SCS basic represents a new standard among laboratory gas burners. With its innovative design this laboratory gas burner is suitable for all flame related applications in the laboratory.

Modern Design

The unique shape is eye-catching. The low profile body facilitates ergonomic operation; the streamlined design reduces air flow disturbances in a clean room workbench to a minimum. Moreover, the Fuego SCS basic is very easy to clean – a laboratory gas burner without grooves and rough edges. Clogging is not an issue. Spilled liquids simply escape through an integrated drain while solid substances can easily be removed by opening the burner shaft. It is even possible to dismantle the individual components of the removable burner head for in-depth cleaning.

Top Rate Safety

The Safety Control System SCS demonstrates state-of-the-art safety technology which constantly analyzes potential hazards and, if necessary, initiates safety measures, such as an interruption of the gas supply. In addition to the ignition and flame control function and the overheating protection, the new safety concept also features continuous burner head control (BHC). BHC unfailingly identifies burner head clogging by liquids or solid

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Fuego SCS basic
Safety laboratory gas burner
with foot pedal

Art.-Nr./Art.-No.
8.001.000

substances and ensures the correct assembly of the burner head. Exceptional passive safety features: After prolonged use, the residual heat display will signal that the burner head is still hot in order to protect the user from burns.

Simple and Safe Operation

The Fuego SCS basic can be operated with a foot pedal. Two programs are available: a flexible start-stop function and the conventional foot pedal operation for short flame sterilization.

Robust

Fabricated entirely of stainless steel, with fireproof controls and a display, protected by heat resistant glass, the Fuego SCS basic can withstand extreme lab conditions.

Fuego SCS basic - Ideal for use in cleanroom workbenches and laboratories.

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The range:

Fuego SCS basic Art.-No. 8.001.000

with foot pedal
2 standard-programs for foot pedal
SCS (Safety Control System) with BHC (Burner Head Control)
Removable and decomposable burner head
Tilt mechanism, right / left
Holding device for 3 inoculation loop holders
Nozzles for natural gas, propane / butane gas
Turbo flame
Wrench 17 mm for gas connection
Screwdriver for burner head and cover of burner shaft
Tubing connector with swivel nut
Switching power supply (global)
Instruction manual and 2-year warranty

Technical Data:

Technology	Microprocessor
Programs	
Foot pedal:	Standard (flame during pressed foot pedal) Start-Stop with timer, 60 min

Safety features

Safety Control System (SCS)
with gas safety cut off: ignition and flame control
temperature monitor
burner head clogging
and assembly monitor (BHC)
automatic unit switch off, 4 hrs.
residual heat display

Gas supply and consumption

Gas connection:	1/4" left with gas filter
Gas types:	II2ELL3B/P: natural gas E/LL, 18 - 25 mbar liquid gas, 20 - 50 mbar 70 g/h liquid gas
Connected load:	CV 360 - 40 min, Express 444 - 50 min, CG 1750 - 150 min, C 206 - 170 min, CP 250 - 210 min, CV 470 - 370 min

Temperatures

Flame temperature:	1350 °C on liquid gas 1300 °C on natural gas (E)
Temperature threshold level:	1.2 kW liquid gas, 1kW natural gas

Electrical

Power consumption:	2 VA
Power connection:	100 - 240 V / 50/60 Hz / max. 0.3 A 9 V DC / 1 A

Mechanical

Casing and operating controls:	stainless steel / glass, UV and solvent resistant
Burner head:	removable and decomposable, stainless steel Ø 23 mm, with drains
Cover of burner shaft:	Measurements (w x h x d): 103 x 49 x 130 mm
Weight:	700 g

Licenses

DIN-DVGW Reg.-No.:	NG-2211AS0167
CE:	EN 61326-1, EN61000-3-2, EN 61010
EEC guidelines:	89/336/EEC, 73/23/EEC

Turboflamme
im Handumdrehen eine extra
harte und stabile Flamme
Turbo flame
an extremely firm and
consistent flame with only a few
flicks of the wrist

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WLD-TEC GmbH