



Home > Stains > Haematology Stains > Special buffers for haematology



Special buffers for haematology

Supplier: VWR Chemicals

[\(No Reviews\)](#)



Buffer tablets Gurr® for the preparation of microscopical stain solutions. These tablets provide a rapid and convenient method of preparing buffer solution. One tablet dissolved in distilled water and made up to 100 ml produce a solution of the stated pH at 20 °C. The molarity of the buffer solutions is 0,005. Contains potassium, sodium and phoshate.

RELATED INFORMATION

Certificates

Product specification

spec_VWRC331552S

spec_VWRC331942F

spec_VWRS331553T

spec_VWRC331932D

SDS

Click to enlarge

