


**Standard Form Crucibles**

CODE NUMBER	CAPACITY CC	DIAMETER mm	DEPTH mm	ESTIMATED WEIGHT (g)					DIAGRAM
				STANDARD FORM	RIM ONLY	BASE ONLY	RIM & BASE	STANDARD LID ONLY	
GCS2.5	2.5	12	15	3	4	4	5	2	
GCS5	5	20	22	4	5	5	6	2	
GCS8	8	25	24	6	7	7	8	2	
GCS10	10	25	28	7	7.5	8	9	3	
GCS15	15	30	31	11	12	13	14	4	
GCS20	20	33	35	15	16	17	18	5	
GCS25	25	36	34	19	20	21	23	5	
GCS30	30	40	36	22	24	26	27	7	
GCS35	35	41	38	25	26	27	32	7	
GCS40	40	42	42	30	32	34	36	7	
GCS50	50	45	44	38	40	42	45	8	
GCS60	60	45	47	46	48	50	54	9	
GSC70	70	49	53	50	52	54	63	12	
GCS80	80	51	53	51	54	63	72	14	
GCS90	90	54	57	52	59	71	81	15	
GCS100	100	56	57	59	68	73	90	17	
GCS110	110	56	60	68	71	81	95	19	
GCS120	120	64	66	84	88	90	98	22	
GCS150	150	65	69	100	104	107	114	27	
GCS200	200	80	65	150	159	165	180	30	
GCS250	250	81	73	175	179	183	190	32	
GCS700	700	107	105	300	305	315	320	50	



**Available Capacity:**  
 2.5-700ml (larger sizes made to order)  
**Options:**  
 Lids, reinforced rims and/or bases.