

## BagPipet® & BagTips®

### Pipet and sterile tips

BagPipet® is an autoclavable pipet with fixed volume (0.1 mL / 0.9 mL / 1 mL). It is used with BagTips®, a single-use sterile and rigid straw.

*BagTips®: its length reaches the bottom of the bag*

*Autoclavable and calibrated according to the GLP*



#### Advantages

- ▶ Economical
- ▶ No risk of cross-contamination
- ▶ Variable volume (0.1 mL/0.9 mL/1 mL)
- ▶ Light and ergonomic
- ▶ Autoclavable

BagTips® is available in 2 sizes:  
Regular (19 cm)/Jumbo (24 cm)


#### Technical specifications

	BagPipet®	BagTips® (Regular)	BagTips® (Jumbo)
Size		190 mm	240 mm
Reference	251 091	252 019	252 024

## BagRack®

### Storage rack for bags

BagRack® enables clean and methodical storage of bags.



400 mL  
Security driptray in option

#### Advantages

- ▶ Adapted to all bags
- ▶ Stackable
- ▶ Robust
- ▶ Stainless steel
- ▶ Autoclavable

Available for bag sizes:  
80 mL, 100 mL, 400 mL, 3500 mL

#### Technical specifications

BagRack®	80/100 mL	400 mL	3500 mL 4 compartments	3500 mL 2 compartments
Size (w x d x h)	12.5 x 27.5 x 14 cm	20.5 x 35.5 x 20 cm	43 x 35 x 37.5 cm	35.5 x 20.5 x 20 cm
Reference	221 010	221 040	221 350	221 352

BagClip®	80/100 mL	400 mL	3500 mL
Quantity per pack	50	50	50
Reference	231 010	231 040	231 350

## BagClip®

### Closing clip

BagClip® allows a temporarily closing in 3 points of any type of bag.



#### Advantages

- ▶ Handy and easy-to-use
- ▶ Gas and liquid leakproof
- ▶ No contact with the sample
- ▶ Avoids intermediary weldings
- ▶ Robust and long-lasting