



#### **Speed control and precision**

Intuitive and convenient speed adjustment simply with the tips of your fingers. The new valve technology enables more sensitive control of the meniscus resulting in increased precision.



#### **Ergonomics**

Lightweight, well-balanced and ergonomically improved design allows for stress-free work. Following our PhysioCare Concept®, we strive to improve your time at work.



#### **Battery status display**

Vibrant backlit LEDs provide optical feedback of the remaining battery life every time the Easypet 3 is used, giving you this extra bit of information just when you need it.

#### Technical Data

Rechargeable battery type	Lithium-polymer
Capacity	1100 mAh/3.7 V
Charging time	~3 h
Number of dispensings	~2000 (with a 25-mL pipette)
Weight without rechargeable battery and pipette	134 g (0.295 lb)
<b>Power supply</b>	
Input	100 V – 240 V AC, $\pm 10\%$ , 50 Hz – 60 Hz, 0.5 A
Output	5 V DC, 1.0 A

#### Ordering Information

Description	International Order No.	Order No. North America
Eppendorf Easypet® 3 incl. power supply, wall mounting device, shelf stand, 2 membrane filters (unsterile) 0.45 $\mu\text{m}$	4430 000.018	–
Eppendorf Easypet® 3 incl. power supply, wall mounting device, 2 membrane filters (unsterile) 0.45 $\mu\text{m}$	–	4430000026
Membrane filter, sterile, 0.45 $\mu\text{m}$ , pack of 5	4421 601.009	022232002
Membrane filter, sterile, 0.2 $\mu\text{m}$ , pack of 5	4430 606.005	4430606005
Lithium-polymer rechargeable battery for Eppendorf Easypet® 3	4430 605.009	4430605009

Your local distributor: [www.eppendorf.com/worldwide](http://www.eppendorf.com/worldwide)

Eppendorf AG · Hamburg · Germany · Tel: +49 40 538 01-0

Application Support E-mail: [support@eppendorf.com](mailto:support@eppendorf.com)

Eppendorf North America, Inc. · USA · Tel: +1 516 334 7500

Toll free phone: +1 800 645 3050 · [techserv@eppendorf.com](mailto:techserv@eppendorf.com)

[www.eppendorf.com/easypet](http://www.eppendorf.com/easypet) · [www.eppendorf.com/tcc](http://www.eppendorf.com/tcc)

eppendorf®, Easypet® and PhysioCare Concept® are registered trademarks of Eppendorf AG, Hamburg, Germany.  
All rights reserved, including graphics and images. Copyright © 2012 by Eppendorf AG.  
Order No.: A442039021/GBZ/3.5T/0610/FEEL/STEF