



The Slide Drying Bench aids in the preparation of microscope slides at the specimen mounting stage. It complements both the Paraffin Wax Dispenser and the Paraffin Wax Mounting Bath to provide a complete set of paraffin was sampling equipment for the modern laboratory.

Accepting up to 50 slides, it has the facility for drying slides in different orientations- across the drying support bars, angled from the bars, or flat without the bars.

There is an On/Off switch and "Power On" white neon and "Heater On" amber neon indicators.

The silicone rubber mat element heater provides heating up to 100°C, and the temperature is controlled from 10°C to 100°C by a built-in energy regulator.

The case and top are powder-coated aluminium and the unit comes coplete with handy carrying handles.

### Features

- Up to 50 slides capacity
- The slides can be arranged on the bench in different orientations
- Built-in energy regulator controlling temperature up to 100°C
- 2 Neon indicators: Clear to indicate "Power On" and amber to indicate "Heater On"
- Replacement elements easy fitted by user

Warranty: One year parts and labour

Certification: CE

### Technical Specification

Material	Base and top powder- coated aluminium
Max Element Temperature	100°C
Dimensions (d x w x h), mm	180 x 390 x 95
Weight, kg	1.8

### Ordering Information

Model	Capacity	Electrical Requirement
MH6616	50 slides	230V, 50/ 60Hz, 150W