

BagMixer® 400 400 mL lab blenders

World best-seller in labs, BagMixer® 400 is an efficient lab blender for all microbiological analyses.



BagMixer® 400 P

- Fixed speed (8 strokes/sec.)
- Variable time (30 - 210 secs. or ∞)
- Also available in CoolMix® version



BagMixer® 400 W

- Fixed speed (8 strokes/sec.)
- Variable time (30 - 210 secs. or ∞)
- **Window door**
- **Security drip tray**

Best seller



BagMixer® 400 VW

- Window door
- Security drip tray
- **Multifunction digital screen**
- **Variable speed (4, 6, 8, 10 strokes/sec.)**
- **Variable time (1 sec. - 1 h or ∞/count-down)**
- **Variable power, adjustable paddles (patented)**
- **Q-Tight® closing**
- **Side-by-side paddle stop**
- **Silent brushless motor**



BagMixer® 400 CC®

- Window door
- Security drip tray
- Multifunction digital screen
- Variable speed (4, 6, 8, 10 strokes/sec.)
- Variable time (1 sec. - 1 h or ∞/count-down)
- Variable power, adjustable paddles (patented)
- Q-Tight® closing
- Side-by-side paddle stop
- Silent brushless motor
- **Click & Clean® removable paddles (patented)**

Demo video
www.interscience.fr

3
YEAR
GUARANTEE

WINDOW DOOR
LIFE-TIME
GUARANTEE

ISO
7218

SHOCK ABSORBERS
LIFE-TIME
GUARANTEE

Lab blenders

Sample preparation

Applications

- Homogenization of all types of samples up to 40 g
- Microbiology: food, pharmaceutical, clinical, veterinary, environmental...

Advantages

- ▲ Silent and robust
- ▲ Variable blending power
- ▲ Easy cleaning: door opens at 270°

Discover the BagMixer® accessories p. 20!

Delivered with: Power cord, user's manual, security drip tray*, free sterile packs of filter bags. (*in option for BagMixer® 400 P)

Technical specifications

BagMixer® 400	P	W	VW	CC®
Useful volume	50 - 400 mL	50 - 400 mL	50 - 400 mL	50 - 400 mL
Power	110/220V/- 50-60Hz	110/220V/- 50-60Hz	100-240V/- 50-60Hz	100-240V/- 50-60Hz
Size (w x d x h)	26 x 39 x 29 cm	26 x 39 x 29 cm	26.5 x 42 x 26 cm	26.5 x 42 x 26 cm
Weight	16.5 kg	16.5 kg	16.5 kg	16 kg
Reference	021 230	022 230	023 230	024 230

Ideal with the sterile bags

