

Messkolben aus DURAN, Klasse DIN A, Ringmarke und Beschriftung *coloraür* blau, mit Bördelrand


Volumetric Flasks DURAN, Class DIN A, graduations *coloraür* blue, with bondering edge



Größe Size	Klasse DIN A	Duran Klarglas Bördelrand	Duran Braunglas Bördelrand
	Class DIN A	Duran clearglass Bondering edge	Duran AMBER Bondering edge
	Toleranz Tolerance	Best.-Nr.: Ord.-No.:	Best.-Nr.: Ord.-No.:
5 ml	± 0,025 ml	0.630.01	0.631.01
10 ml	± 0,025 ml	.02	.02
20 ml	± 0,040 ml	.03	.03
25 ml	± 0,040 ml	.04	.04
50 ml	± 0,060 ml	.05	.05
100 ml	± 0,100 ml	.06	.06
200 ml	± 0,150 ml	.07	.07
250 ml	± 0,150 ml	.08	.08
500 ml	± 0,250 ml	.09	.09
1000 ml	± 0,400 ml	.10	.10
2000 ml	± 0,600 ml	.11	.11

Messkolben aus DURAN, Klasse DIN A, Ringmarke und Beschriftung *coloraür* blau, mit NS-K-Stopfen

Volumetric Flasks DURAN, Class DIN A, graduations *coloraür* blue, with NS-plastic stopper



Größe Size	Klasse DIN A	Duran Klarglas mit NS K-Stopfen	Duran Braunglas mit NS K-Stopfen
	Class DIN A	Duran clearglass NS-plastic stopper	Duran AMBER NS-plastic stopper
	Toleranz Tolerance	Best.-Nr.: Ord.-No.:	Best.-Nr.: Ord.-No.:
1 ml NS 7	± 0,025 ml	DURAN	0.632.01
2 ml NS 7	± 0,025 ml	DURAN	.02
5 ml NS 10	± 0,025 ml	DURAN	.03
10 ml NS 10	± 0,025 ml	DURAN	.04
20 ml NS 10	± 0,040 ml	DURAN	.05
25 ml NS 10	± 0,040 ml	DURAN	.06
50 ml NS 12	± 0,060 ml	DURAN	.07
100 ml NS 12	± 0,100 ml	DURAN	.08
200 ml NS 14	± 0,150 ml	DURAN	.10
250 ml NS 14	± 0,150 ml	DURAN	.11
500 ml NS 19	± 0,250 ml	Boro 3.3	.12
1000 ml NS 24	± 0,400 ml	Boro 3.3	.13
2000 ml NS 29	± 0,600 ml	Boro 3.3	.14
5000 ml NS 34	± 1,200 ml	DURAN	.15
10000 ml NS 45	± 2,400 ml	DURAN	.16
500 ml NS 19	± 0,250 ml	DURAN	.20
1000 ml NS 24	± 0,400 ml	DURAN	.21
2000 ml NS 29	± 0,600 ml	DURAN	.22