


Titrierbüretten nach Pellet Klasse DIN AS, ohne Zwischenhahn, Duran, Klarglas, Hauptpunkte Ringteilung *colordur* blau, ohne Flasche


Titration Burettes Pellet Class DIN AS, without stopcock between burette and bottle, Duran-clearglass, main points ring graduation *colordur* blue, without bottle




Größe Size	Klasse DIN AS	Duran, Klarglas NS-Glashahn	Duran, Klarglas NS-PTFE Hahn	Duran, Klarglas PTFE- Spindelhahn
	Class DIN AS	Duran, clearglass NS-glass stopcock	Duran, clearglass NS-PTFE stopcock	Duran, clearglass PTFE-spindle valve
	Toleranz Tolerance	Best.-Nr.: Ord.-No.:	Best.-Nr.: Ord.-No.:	Best.-Nr.: Ord.-No.:
10 ml : 1/50	± 0,02 ml	0.552.01	0.552.11	0.552.21
10 ml : 1/20	± 0,02 ml	.02	.12	.22
25 ml : 1/20	± 0,05 ml	.03	.13	.23
25 ml : 1/10	± 0,05 ml	.04	.14	.24
50 ml : 1/10	± 0,10 ml	.05	.15	.25
100 ml : 1/5	± 0,20 ml	.06	.16	.26

Titrierbüretten nach Pellet Klasse DIN AS, ohne Zwischenhahn, Duran, Schellbach, Hauptpunkte Ringteilung *colordur* blau, ohne Flasche

Titration Burettes Pellet Class DIN AS, without stopcock between burette and bottle, Duran-Schellbach, main points ring graduation *colordur* blue, without bottle



Größe Size	Klasse DIN AS	Duran, Schellbach NS-Glashahn	Duran, Schellbach NS-PTFE Hahn	Duran, Schellbach PTFE- Spindelhahn
	Class DIN AS	Duran, Schellbach NS-glass stopcock	Duran, Schellbach NS-PTFE stopcock	Duran, Schellbach PTFE-spindle valve
	Toleranz Tolerance	Best.-Nr.: Ord.-No.:	Best.-Nr.: Ord.-No.:	Best.-Nr.: Ord.-No.:
10 ml : 1/50	± 0,02 ml	0.553.01	0.553.11	0.553.21
10 ml : 1/20	± 0,02 ml	.02	.12	.22
25 ml : 1/20	± 0,05 ml	.03	.13	.23
25 ml : 1/10	± 0,05 ml	.04	.14	.24
50 ml : 1/10	± 0,10 ml	.05	.15	.25
100 ml : 1/5	± 0,20 ml	.06	.16	.26



Zubehör für GLASFIRN Titrierbüretten nach Pellet

Accessories to GLASFIRN Titration Burettes Pellet

Best.-Nr.:

Ord.-No.:

0.560.10	Bürettenflasche Klarglas 1 L NS 29
.20	Bürettenflasche Klarglas 2 L NS 29
0.562.10	Bürettenflasche Braunglas 1 L NS 29
.20	Bürettenflasche Braunglas 2 L NS 29
0.564.10	Gummidruckball

0.560.10	Reservoir Flasks clear 1 L NS 29
.20	Reservoir Flasks clear 2 L NS 29
0.562.10	Reservoir Flasks amber 1 L NS 29
.20	Reservoir Flasks amber 2 L NS 29
0.564.10	Rubber bulb