

AUTOCLAVE BASKET

Description	Square basket
Dimension, W x D x H, mm, approx.	200 x 200 x 200
Size of grid / mm, approx.	3 x 3
Material	Stainless steel type 304
Application	To store and carry materials for safe autoclaving
Product code	BG230-30

Description	Square basket
Dimension, W x D x H, mm, approx.	150 x 150 x 150
Size of grid, mm, approx.	3 x 3
Material	Stainless steel type 304
Application	To store and carry materials for safe autoclaving
Product code	BG230-15

Description	Square basket
Dimension, W x D x H, mm, approx.	100 x 100 x 100
Size of grid, mm, approx.	2 x 2
Material	Stainless steel type 304
Application	To store and carry materials for safe autoclaving
Product code	BG230-10

