

DRYING CABINET

Description	Complete with 2 perforated height adjustable shelves, 2 sliding glass doors, SIRIM-approved safety thermostatic control, and 13A safety fused plug top. Up to 80°C.
Dimension, W x D x H, mm	860 x 500 x 1000
Material of the Cabinet	Light grey epoxy-coated mild steel, single wall
Material of the Shelf	Stainless steel type 304
Application	To rapidly dry glasswares or samples.
Product code	OV100-85

Note:

Custom-made Cabinet is available upon request.

