

HOTPLATE MAGNETIC STIRRER (manufactured from South Korea)

Description	Top plate is resistant to most chemicals and thermal crack. Comes with stainless steel lattice rod. 1 knob to control stirring speed and 1 knob to control the heating temperature.
Stirring speed, rpm	60 to 1500
Heating temperature, °C	60 to 420
Stirring capacity	2 litre
Exterior Dimensions, W x D x H, mm Dimensions of the plate, W x D, mm	205 x 260 x 110 190 x 190
Material of top plate	Ceramic coated, stainless steel
Weight, kg, approximate	3.5
Application	To heat and stir solution in a vessel, all at once.
Warranty	1 year against manufacturing defect
Product code	HTS-1003

Note :

Please refer to brand “Assistent, Germany” if you wish to know more about magnetic stir bar and retriever.

Please refer to brand Kimble Chase and GQ if you wish to know more about glasswares.

