


PETRI DISH CARRIER

Description	<p>One wired compartment</p> <p>Comes with a handle</p> <p>Suitable for glass and plastic Petri dish</p>
Dimension, I.D. x H, mm, approx	102 x 220
Material	Epoxy-coated mild steel
Application :	<p>To carry Petri dish for safe autoclaving with the following quantity:</p> <p>15 x Petri dish (90 mm dia. x 15 mm H)</p> <p>10 x Petri dish (100 mm dia. x 20 mm H)</p>
Product code	PD410-10



Labmart

Description	3 solid compartments Comes with a handle Suitable for glass and plastic Petri dish
Dimension, I.D. x H, mm, approx	102 x 220 per compartment
Material	Epoxy-coated mild steel
Application :	To carry Petri dish for safe autoclaving with the following quantity: 15 x Petri dish (90 mm dia. x 15 mm H) per compartment x 3 10 x Petri dish (100 mm dia. x 20 mm H) per compartment x 3
Product code	PD430-10

Note :

Autoclavable (121 °C, 15 psi/ 100 kPa above atmospheric pressure).