



CORROSIVE CHEMICAL SAFETY STORAGE CABINET (Blue)

| | |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Complete with 1 perforated fixed height tray, 1 stainless steel spill collection tray, 2 swing-out doors with lock, safety signage in dual languages, and mounted on 4x4" nylon wheels with brake. * Additional tray can be purchased separately. |
| Dimension, W x D x H, mm | 750 x 600 x 820 |
| Material of the Cabinet | Blue epoxy-coated mild steel |
| Application | To safely and professionally store various corrosive chemicals (for example, acids, alkaline, etc) |
| Product code | SBB7530-A |





| | |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Complete with 3 height adjustable trays, 1 stainless steel spill collection tray, 2 swing-out doors with lock, safety signage in dual languages, and mounted on 4x4" nylon wheels with brake. * Additional tray can be purchased separately. |
| Dimension, W x D x H, mm | 850 x 600 x 1625 |
| Material of the Cabinet | Blue epoxy-coated mild steel |
| Application | To safely and professionally store various corrosive chemicals (for example, acids, alkaline, etc) |
| Product code | SBB1625-A |

Note :

Custom-size Cabinet is available upon request, including stainless steel material.