



CHEMICAL SAFETY STORAGE CABINET (Yellow)

Description	Complete with 1 perforated fixed height shelf, 1 stainless steel spill collection tray, 2 swing-out doors with lock, safety signage in dual languages, and mounted on 4x4" nylon wheels with brake. * Additional shelf can be purchased separately.
Dimension, W x D x H, mm	750 x 600 x 820
Material of the Cabinet	Yellow epoxy-coated mild steel
Application	To safely and professionally store various flammable chemicals (for example, ethanol, acetone, xylene etc)
Product code	SBB7530





Description	<p>Complete with 3 perforated height adjustable shelves, 1 stainless steel spill collection tray, 2 swing-out doors with lock, safety signage in dual languages, and mounted on 4x4" nylon wheels with brake.</p> <p>* Additional shelf can be purchased separately.</p>
Dimension, W x D x H, mm	850 x 600 x 1625
Material of the Cabinet	Yellow epoxy-coated mild steel
Application	To safely and professionally store various flammable chemicals (for example, ethanol, acetone, xylene etc)
Product code	SBB1625

Note :

Custom-size Cabinet is available upon request, including stainless steel material.