

SLIDE STORAGE CABINET

Description	<p>Comes with 14 vertical pull-out drawers</p> <p>Each drawer comes a cylindrical chromed steel handle</p> <p>Compact and space-saving design</p>
Overall Dimension, W x D x H, mm (approximate)	480 x 480 x 140
Dimension of Drawer, W x D x H, mm (approximate)	30 x 460 x 120
Material	Light grey epoxy-coated mild steel
Application	To neatly organise and store microscopic slides 3"x1"
Product code	SC14F



Description	<p>Comes with 10 horizontal pull-out drawers</p> <p>Each drawer comes with ready-slots for slides, an intact frame for labelling, and a chromed steel handle.</p> <p>Compact and space-saving design</p>
Overall Dimension, W x D x H, mm (approximate)	360 x 550 x 420
Dimension of Drawer, W x D x H, mm (approximate)	340 x 520 x 38
Material	Beige epoxy-coated mild steel
Application	To neatly organise and store microscopic slides 3"x1"
Product code	SC10D