

STERILISATION CANISTER – GLASS PIPETTE

Description	Cylindrical container used for Pasteur pipette Suitable for Pasteur pipette : 150 mm and 230 mm
Dimension, Height x Diameter, mm, approx	290 x 60
Material	Stainless steel type 304
Application	To store Pasteur / transfer pipette for safe autoclaving.
Product code	PP970-15

Description	Cylindrical container used for glass pipette Suitable for : Measuring pipette, 1-20 ml Volumetric / Bulb pipette, 1-7 ml
Dimension, Height x Diameter, mm	400 x 60
Material	Stainless steel type 304
Application	To store measuring pipette & bulb/volumetric pipette for safe autoclaving.
Product code	PP970-35



Description	Cylindrical container used for glass pipette Suitable for Measuring pipette, 25 ml and 50 ml; Volumetric / Bulb pipette, 15 ml
Dimension, W x D x H, mm	500 x 60
Material	Stainless steel type 304
Application	To store measuring pipette & bulb/volumetric pipette for safe autoclaving.
Product code	PP970-55

Note :

Autoclavable (121 °C, 15 psi/ 100 kPa above atmospheric pressure).