

LABORATORY TROLLEY

| | |
|--------------------------|--|
| Description | <p>2 tiers, with 1" high flip up sides, reinforced lower tier support</p> <p>Max weight capacity : 75 kg / tier</p> <p>Complete with 4 x 4" rubber castor wheels with brake.</p> <p>Rust resistance.</p> |
| Dimension, W x D x H, mm | 700 x 450 x 900 |
| Material | Stainless steel type 304 |
| Application | To conveniently cart and transport various laboratory apparatus and consumables. |
| Product code | TT800-12 |

| | |
|--------------------------|--|
| Description | <p>3 tiers, with 1" high flip up sides, reinforced lower tier support</p> <p>Max weight capacity : 75 kg / tier</p> <p>Complete with 4 x 4" rubber castor wheels with brake.</p> <p>Rust resistance.</p> |
| Dimension, W x D x H, mm | 700 x 450 x 900 |
| Material | Stainless steel type 304 |
| Application | To conveniently cart and transport various laboratory apparatus and consumables. |
| Product code | TT800-16 |

Labmart

LABORATORY TROLLEY

