

## VORTEX MIXER

Description	<p>Simple operation with ON/OFF switch and a rotary knob for speed adjustment.</p> <p>Comes with 2 adapters: 1 adapter for tubes (up to 50ml), and the other adapter for flat bottom flask (up to 500ml).</p> <p>4 studs at the bottom all corners to provide steady footing.</p>
Adjustable Speed	20 to 3000 rpm
Exterior Dimensions, W x D x H, mm	125 x 150 x 165
Material	ABS body
Application	To mix liquid in small vessel, or to provide circular agitation to the liquid in the vessel
Product code	LM-3000 (formerly coded as VTX-3000L) (fully made in South Korea)

