


Büretten mit geradem Hahn, Klasse DIN AS, DURAN, Klarglas oder Klarglas mit Schellbachstreifen Hauptpunkte Ringteilung **colordur** blau


Burettes with straight stopcock, Class DIN AS, clearglass or clearglass with Schellbach-strip, main points ring-Graduation **colordur** blue



	Klasse DIN AS	Duran, Klarglas NS-Glashahn	Duran, Schellbach NS-Glashahn	Duran, Klarglas NS-PTFE Hahn	Duran, Schellbach NS-PTFE Hahn
	Class DIN AS	Duran, clearglass NS-glass stopcock	Duran, Schellbach NS-glass stopcock	Duran, clearglass NS-PTFE stopcock	Duran, Schellbach NS-PTFE stopcock
Größe Size	Toleranz Tolerance	Best.-Nr.: Ord.-No.:	Best.-Nr.: Ord.-No.:	Best.-Nr.: Ord.-No.:	Best.-Nr.: Ord.-No.:
5 ml : 1/50	± 0,01 ml	0.582.01	0.583.01	0.582.11	0.583.11
10 ml : 1/50	± 0,02 ml	.02	.02	.12	.12
10 ml : 1/20	± 0,02 ml	.03	.03	.13	.13
10 ml : 1/10	± 0,02 ml	.04	.04	.14	.14
25 ml : 1/10	± 0,05 ml	.05	.05	.15	.15
25 ml : 1/20	± 0,05 ml	.06	.06	.16	.16
50 ml : 1/10	± 0,10 ml	.07	.07	.17	.17
100 ml : 1/5	± 0,20 ml	.08	.08	.18	.18

Büretten mit geradem Hahn, Klasse DIN AS, DURAN-BRAUNGLAS, Hauptpunkte Ringteilung **colordur** weiß

Burettes with straight stopcock, Class DIN AS, DURAN-AMBER-GLASS, main points ring-graduation **colordur** white



	Klasse DIN AS	Duran, BRAUNGLAS NS-Glashahn	Duran, BRAUNGLAS NS-PTFE Hahn
	Class DIN AS	Duran, AMBER NS-glass stopcock	Duran, AMBER NS-PTFE stopcock
Größe Size	Toleranz Tolerance	Best.-Nr.: Ord.-No.:	Best.-Nr.: Ord.-No.:
5 ml : 1/50	± 0,01 ml	0.584.01	0.584.11
10 ml : 1/50	± 0,02 ml	.02	.12
10 ml : 1/20	± 0,02 ml	.03	.13
25 ml : 1/20	± 0,05 ml	.05	.15
50 ml : 1/10	± 0,10 ml	.06	.16
100 ml : 1/5	± 0,20 ml	.07	.17