

RESTRAINER

Description	<p>Mice restrainer</p> <p>Comes with 3 slots for different sizes (length) of mice.</p>
Exterior Dimensions, W x L x H, mm, approx	40 x 140 x 40
Material	Clear perspex / acrylic
Application	To restrain mice from moving about so that further action can be taken on the mice, for instance, injection.
Product code	R5140



Slots


Description	<p>Rat restrainer</p> <p>Comes with 4 slots for different sizes (length) of rats.</p>
Exterior Dimensions, W x L x H, mm, approx	60 x 200 x 60
Material	Clear perspex / acrylic
Application	To restrain rat from moving about so that further action can be taken on the rat, for instance, injection .
Product code	R5200



RESTRAINER

Description	Rabbit restrainer Comes with 3 holes to fasten the lock to adjust the holes for the rabbit neck size.
Exterior Dimensions, W x L x H, mm, approx	160 x 400 x 220
Material	Clear perspex / acrylic
Application	To restrain rabbit from moving about so that further action can be taken on the rabbit, for instance, injection.
Product code	R5300

Lock



Hole to restrain the
rabbit neck