



## Product Specification Sheet

# VWR SPILL RESPONSE KIT FOR CHEMICALS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

*VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Cat. No.: **121-1050**

Dimensions: **610 x 150 x 480mm (w x d x h)**

Product weight: **1.6Kg**

Absorbency: up to 10 L, depending on the chemical being absorbed

Use to absorb: **Oils, fuels, solvents, cutting fluids, aggressive acids, aggressive alkalis, milk, water, etc.,**

### Kit contents:

|                                                            |     |
|------------------------------------------------------------|-----|
| 7.5cm x 1.2m 'Superior' maintenance sock, capacity 3.75l   | 1 x |
| 51cm x 48cm maintenance Absorbent pad, capacity 0.9 l/each | 9 x |
| 30cm black cable tie                                       | 2 x |
| 46cm x 90cm 200 gauge blue plastic disposal bag            | 2 x |
| 'Silver Shield' Chemical resistant glove                   | 1 x |
| 60cm x 48cm clear plastic clip-closed carrier              | 1 x |
| Non branded spill kit instruction sheet                    | 1 x |

**Disclaimer:** VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt ask us.

Date of Issue: **August, 2011**